

FAHLL

AIRWAY  
MANAGEMENT



MEDICAL DEVICES FOR  
LARYNGECTOMY AND  
TRACHEOTOMY CARE

**Edition/date: 03/2019**

The manufacturer reserves the right to alterations to the products listed. especially technical alterations.

### Copyright

The contents of this catalogue. including all of its parts and images are protected by copyright. Any use beyond the narrow limits of the copyright without the consent of the author is prohibited and subject to prosecution. This applies. in particular. to reproductions. translations. microfilming and the storage and processing of the data in electronic systems.

This brochure contains only a small selection of our product range. More about our range of products you will find at [www.fahl.de](http://www.fahl.de).

## LOCATIONS

### Andreas Fahl Medizintechnik-Vertrieb GmbH

✉ August-Horch-Str. 4a  
51149 Köln - Germany  
☎ +49 (0) 22 03/29 80-0  
☎ +49 (0) 22 03/29 80-100  
✉ info@fahl.de  
🌐 www.fahl.de

### Local branch Berlin

✉ Am Borsigturm 62  
13507 Berlin - Germany  
☎ +49 (0) 30/43 73 139-0  
☎ +49 (0) 30/43 73 139-35  
✉ info@fahl.de  
🌐 www.fahl.de

### Overseas branch Austria

✉ Triester Straße 14 / Bürohaus 2  
2351 Wiener Neudorf - Austria  
☎ +43 (0) 22 36/89 28 38-0  
☎ +43 (0) 22 36/89 28 38-14  
✉ info@fahl.de  
🌐 www.fahl.de

### Fahl Scandinavia AB

✉ Building 301 Medicon Village  
223 81 Lund - Schweden  
☎ +46 (0) 46 - 287 37 90  
✉ info@fahl-scandinavia.se  
🌐 www.fahl-scandinavia.se


### Local branch Munich





















✉ Aurbacher Straße 2  
81541 München - Germany  
☎ +49 (0) 89/58 90 90 9-0  
☎ +49 (0) 89/58 90 90 9-20  
✉ info@fahl.de  
🌐 www.fahl.de

### Overseas branch Switzerland

✉ Cherstrasse 4  
8152 Glattbrugg - Schweiz  
☎ +41 (0) 448293060  
☎ +41 (0) 448293061  
✉ info@fahl-medical.ch  
🌐 www.fahl-medical.ch

## HOW TO FIND US - HEAD OFFICE COLOGNE



















TRACHEOSTOMY TUBES		4
TRACHEOSTOMY TUBES - LINGO / PHON		14
TRACHEOSTOMY TUBES - CUFF		22
TRACHEOSTOMY TUBES - LINGO CUFF / PHON		27
TRACHEOSTOMY TUBES - SUCTION		33
VOICE PROSTHESES		36
SPEAKING VALVES		41
TRACHEOSTOMA VALVES		43
HME FILTER CASSETTES / BASE PLATES		44
TRACHEAL COMPRESSES		53
CLEANING AND CARE OF CANNULAS		56
TUBE HOLDERS		59
TRACHEOSTOMA PROTECTION		63
INHALATION AND AIR HUMIDIFICATION		68
TRACHEAL SUCTION DEVICES		70
SPEECH AIDS / VOICE AMPLIFIERS		72
EMERGENCY AIDS		73
WATER THERAPY		75
CLINICAL PRODUCTS		76
ANATOMY / FURTHER INFORMATION		79









# TRACHEOSTOMY TUBES

















We provide a wide range of tracheostomy tubes made from various materials and with different design patterns. This diversity is due to the fact that tracheostomy tubes should be ideally matched to the individual requirements of each patient.















The doctor, who individually selects and adjusts the tubes for each patient, should be made aware of the complete range to allow him to select the best product for his patients.

	Description	Unit / pcs	Size	REF	
 	<b>SILVERVENT®</b> 1 IC	1	1 - 13	<b>10000</b>	
 	<b>SILVERVENT® UNI</b> 1 ICU	1	0 - 12	<b>10500</b>	
 	<b>DURAVENT®</b>	1	3 - 13	<b>11000</b>	
 	<b>DURAVENT® XL</b>	1	5 - 13	<b>11200</b>	
 	<b>DURAVENT®</b> 1 IC	1	3 - 13	<b>11001</b>	
 	<b>DURAVENT® XL</b> 1 IC	1	5 - 13	<b>11201</b>	




























	Description	Unit / pcs	Size	REF	
 	<b>DURAVENT®</b> 2 IC	1	3 - 13	<b>11002</b>	
 	<b>DURAVENT® UNI</b> 1 ICU	1	3 - 13	<b>11011</b>	
 	<b>DURAVENT® UNI</b> 1 IC + 1 ICU	1	3 - 13	<b>11022</b>	
 	<b>DURAVENT® UNI XL</b> 1 IC + 1 ICU	1	5 - 13	<b>11222</b>	
 	<b>DURAVENT® UNI</b> 2 ICU	1	3 - 13	<b>11012</b>	
 	<b>DURAVENT® KOMBI</b> 1 IC + 1 ICK	1	5 - 13	<b>11092</b>	
 	<b>DURAVENT® SOLO-ADAPTER KOMBI</b>	5	5 - 13	<b>11990</b>	
 	<b>DURAVENT® SOLO-ADAPTER UNI</b>	5	5 - 13	<b>11992</b>	


	Description	Unit / pcs	Size	REF	
 	DURAPUR®	1	7 - 12	23000	
 	<span style="border: 1px solid black; padding: 2px;">XL</span> DURAPUR® XL	1	8 - 10	23200	
 	DURAPUR® 1 IC	1	7 - 12	23001	
 	<span style="border: 1px solid black; padding: 2px;">XL</span> DURAPUR® XL 1 IC	1	8 - 10	23201	
 	DURAPUR® UNI 1 ICU	1	7 - 12	23011	
 	<span style="border: 1px solid black; padding: 2px;">XL</span> DURAPUR® UNI XL 1 ICU	1	8 - 10	23211	
 	DURAPUR® KOMBI 1 ICK	1	7 - 12	23081	










	Description	Unit / pcs	Size	REF	
 <span data-bbox="443 324 483 365">XL</span> 	<b>DURAPUR® KOMBI XL</b> 1 ICK	1	8 - 10	23281	
 	<b>DURAPUR® IC</b>	1	7 - 12	23801	
 <span data-bbox="443 761 483 801">XL</span> 	<b>DURAPUR® IC XL</b>	1	8 - 10	23701	
 	<b>DURAPUR® ICU</b>	1	7 - 12	23811	
 	<b>DURAPUR® ICK</b>	1	7 - 12	23881	
 	<b>LARYNGOTEC®</b>	1	9 - 12	14120	
 	<b>LARYNGOTEC® KOMBI</b>	1	7 - 13	14920	
 	<b>LARYNGOTEC® KOMBI CLIP</b>	1	7 - 13	14830	











	Description	Unit / pcs	Size	REF	
 	<p><b>DURATWIX® JUNIOR PED UNI</b></p>	1	3 - 6	17401	
 	<p><b>DURATWIX® JUNIOR PED L UNI</b></p>	1	3 - 6	17431	
 	<p><b>SPIRAFLEX® UNI</b></p>	1	6 - 11	17504	
 	<p><b>SPIRAFLEX® MRT UNI</b></p>	1	6 - 11	17521	
 	<p><b>SPIRAFLEX® SHORT UNI 2 ICX</b></p>	1	6 - 11	17506	S
 	<p><b>SPIRAFLEX® ICX</b></p>	3	6 - 11	17543	
 	<p><b>SPIRAFLEX® SHORT ICX</b></p>	3	6 - 11	17516	S







	Description	Unit / pcs	Size	REF	
 <p>Icons:    </p>	<b>TRACHEOTEC® VARIO</b>	1	3 - 10	18560	
 <p>Icons:    </p>	<b>FAHL® BIESALSKI</b> 2 IC + Speaking valve	1	4 - 12	17900	
 <p>Icons:  </p>	<b>FAHL® BIESALSKI IC</b>	3	6 - 12	17905	
 <p>Icons:  </p>	<b>FAHL® BIESALSKI</b> Speaking valve S	2	4 - 7	17901-01	
 <p>Icons:  </p>	<b>FAHL® BIESALSKI</b> Speaking valve M	2	8 - 9	17901-02	
 <p>Icons:  </p>	<b>FAHL® BIESALSKI</b> Speaking valve L	2	10 - 12	17901-03	
 <p>Icons:    </p>	<b>DURATWIX® UNI</b>	1	7 - 10	19201	

	Description	Unit / pcs	Size	REF	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span>    <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® UNI short</b></p>	1	7 - 10	19202	
 <div style="display: flex; align-items: center; gap: 5px;">     <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® UNI-VARIO</b> 1 ICU + 1 ICV</p>	1	6 - 10	19211	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span>     <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® UNI-VARIO short</b> 1 ICU + 1 ICV</p>	1	7 - 10	19212	
 <div style="display: flex; align-items: center; gap: 5px;">  <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® IC</b></p>	3	7 - 10	19801	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span>  <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® IC short</b></p>	3	7 - 10	19802	
 <div style="display: flex; align-items: center; gap: 5px;">  <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® ICU</b></p>	3	6 - 10	19811	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span>  <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® ICU short</b></p>	3	6 - 10	19812	
 <div style="display: flex; align-items: center; gap: 5px;">  <span style="border: 1px solid black; padding: 2px;">STERILE</span> </div>	<p><b>DURATWIX® ICV</b></p>	3	6 - 10	19821	

		Description	Unit / pcs	Size	REF	
	S	DURATWIX® ICV short	3	6 - 10	19822	
		DURATWIX® ICU contrast	3	6 - 10	19803	
		DURATWIX® ICK	3	7 - 10	19831	
	S	DURATWIX® ICK short	3	7 - 10	19832	
		DURATWIX® SOLO-ADAPTER KOMBI	3	22 mm	19808	
		DURATWIX® SOLO-ADAPTER VARIO	3	15 mm	19807	
		DURATWIX® SOLO-ADAPTER UNI	3	15 mm	19806	
		TRACHEOSILC® VARIO	1	5 - 9,5	18820	
		LARYNGOTEC® Stoma Button FIX	1	8 - 12	15900	











		Description	Unit / pcs	Size	REF	
	S	LARYNGOTEC® Stoma Button FIX short	1	8 - 12	15910	
		LARYNGOTEC® Stoma Button	1	8 - 12	15905	
	S	LARYNGOTEC® Stoma Button short	1	8 - 12	15915	
		LARYNGOTEC® KOMBI Stoma Button FIX	1	8 - 12	15950	
	S	LARYNGOTEC® KOMBI Stoma Button FIX short	1	8 - 12	15960	
		LARYNGOTEC® KOMBI Stoma Button	1	8 - 12	15955	
	S	LARYNGOTEC® KOMBI Stoma Button short	1	8 - 12	15965	
		FAHL® T-Tube	1	8 - 14	27502	
		FAHL® Tracheo-Safe C	1	45 mm	20895	
		FAHL® Tracheo-Safe P	1	70 mm	20898	








	Description	Unit / pcs	Size	REF	
	<b>SENSOTRACH® Foam OVAL</b>	10	oval	<b>30790</b>	
	<b>SENSOTRACH® Foam PRO OVAL</b>	10	oval	<b>30791</b>	
	<b>SENSOTRACH® Foam ROUND</b>	10	round	<b>30792</b>	
	<b>SENSOTRACH® Foam PRO ROUND</b>	10	round	<b>30793</b>	










These tracheostomy tubes have an opening in the upper (convex) wall of the tube shaft facing the larynx, the so-called sieve-fenestration. While most conventional tracheostomy tubes have one or two simple holes the Fahl tracheostomy tubes have an elaborate sieve consisting of many small holes to reduce micro-aspiration and the ingrowth of granulation tissue.

If an inner cannula is present, it has an oval opening corresponding to the sieve in the outer cannula. Sieve-fenestrated cannulas are used by tracheostomized patients with intact larynx who are able to speak, or by laryngectomized patients with voice prosthesis who still need a tracheostomy tube to stabilize the trachea and keep the tracheostoma open. In that case the sieve supports the speaking facility provided by the voice prosthesis.

































Sieve-fenestrated tracheostomy tubes are usually used in combination with a speaking valve.

	Description	Unit / pcs	Size	REF	
 	<b>LARYNGOTEC® LINGO</b>	1	9 - 12	<b>14220</b>	
 	<b>LARYNGOTEC® KOMBI LINGO</b>	1	7 - 13	<b>14720</b>	
 	<b>LARYNGOTEC® KOMBI CLIP LINGO</b>	1	7 - 13	<b>14820</b>	
 	<b>LARYNGOTEC® KOMBI CUT</b>	1	7 - 13	<b>14725</b>	
 	<b>PHONTRACH® KOMBI LINGO-PHON</b>	1	7 - 13	<b>14930</b>	

















	Description	Unit / pcs	Size	REF	
 <p>S</p> <p>IK LINGO BE PHON FLEX IK S STERILE</p>	<p><b>SPIRAFLEX® SHORT UNI LINGO-PHON</b> 1 ICX + 1ICFX + Speaking Valve</p>	1	6 - 11	17508	
 <p>S</p> <p>FLEX STERILE</p>	<p><b>SPIRAFLEX® SHORT ICX</b></p>	3	6 - 11	17516	
 <p>S</p> <p>FLEX IK STERILE</p>	<p><b>SPIRAFLEX® SHORT ICFX</b></p>	3	6 - 11	17517	
 <p>STERILE</p> <p>IK LINGO 15mm 15mm IK MRT</p>	<p><b>DURATWIX® UNI-VARIO LINGO</b> 1 ICV + 1 ICFU</p>	1	6 - 10	19261	
 <p>S</p> <p>IK LINGO 15mm 15mm IK MRT</p> <p>STERILE</p>	<p><b>DURATWIX® UNI-VARIO LINGO short</b> 1 ICV + 1 ICFU</p>	1	7 - 10	19262	
 <p>IK LINGO MRT</p> <p>STERILE</p>	<p><b>DURATWIX® UNI LINGO</b></p>	1	7 - 10	19271	
 <p>S</p> <p>IK LINGO MRT</p> <p>STERILE</p>	<p><b>DURATWIX® UNI LINGO short</b></p>	1	7 - 10	19272	



		Description	Unit / pcs	Size	REF	
		DURATWIX® IC	3	7 - 10	19801	
	S	DURATWIX® IC short	3	7 - 10	19802	
		DURATWIX® ICU	3	6 - 10	19811	
	S	DURATWIX® ICU short	3	6 - 10	19812	
		DURATWIX® ICV	3	6 - 10	19821	
	S	DURATWIX® ICV short	3	6 - 10	19822	
		DURATWIX® ICU contrast	3	6 - 10	19803	
		DURATWIX® ICK	3	7 - 10	19831	
	S	DURATWIX® ICK short	3	7 - 10	19832	







	Description	Unit / pcs	Size	REF	
   	DURATWIX® PHON ICF	2	7 - 10	19841	
   	DURATWIX® PHON ICF short	2	7 - 10	19842	S
  	DURATWIX® ICF	3	7 - 10	19851	
  	DURATWIX® ICF short	3	7 - 10	19852	S
  	DURATWIX® ICFU	3	6 - 10	19861	
  	DURATWIX® ICFU short	3	6 - 10	19862	S
  	DURATWIX® ICFV	3	6 - 10	19871	
  	DURATWIX® ICFV short	3	6 - 10	19872	S
  	DURATWIX® ICFK	3	7 - 10	19881	
  	DURATWIX® ICFK short	3	7 - 10	19882	S

	Description	Unit / pcs	Size	REF	
 	<p><b>DURATWIX®</b> <b>SOLO-ADAPTER KOMBI</b></p>	3	22 mm	<b>19808</b>	
 	<p><b>DURATWIX®</b> <b>SOLO-ADAPTER VARIO</b></p>	3	15 mm	<b>19807</b>	
 	<p><b>DURATWIX®</b> <b>SOLO-ADAPTER UNI</b></p>	3	15 mm	<b>19806</b>	
 	<p><b>DURATWIX®</b> <b>PLUG</b></p>	5	15 mm	<b>19805</b>	
    	<p><b>SILVERVENT®</b> <b>LINGO-PHON</b> 1 ICF</p>	1	0 - 14	<b>20000</b>	
      	<p><b>SILVERVENT®</b> <b>MULTI</b> 1 ICU + 1 ICF-O<sub>2</sub></p>	1	2 - 13	<b>20300</b>	
      	<p><b>DURAVENT®</b> <b>UNI-KOMBI</b> <b>LINGO</b> 1 ICFU + 1 ICFK</p>	1	5 - 13	<b>11072</b>	
    	<p><b>DURAVENT®</b> <b>LINGO-PHON</b> 1 IC + 1 ICF</p>	1	3,5 - 13	<b>21032</b>	
     	<p><b>DURAVENT®</b> <b>UNI</b> <b>LINGO-PHON</b> 1 ICU + 1 ICF</p>	1	3,5 - 13	<b>21042</b>	

	Description	Unit / pcs	Size	REF	
 <p><b>XL</b></p> 	<p><b>DURAVENT® UNI LINGO-PHON XL</b> 1 ICU + 1 ICF</p>	1	5 - 13	<b>21242</b>	
 	<p><b>DURAVENT® KOMBI LINGO-PHON</b> 1 IC + 1 ICFK + 1 COMBIPHON®</p>	1	5 - 13	<b>21062</b>	
 	<p><b>DURAVENT® UNI-KOMBI LINGO-PHON</b> 1 ICU + 1 ICFK + 1 COMBIPHON®</p>	1	5 - 13	<b>21072</b>	
 <p><b>XL</b></p> 	<p><b>DURAVENT® UNI-KOMBI LINGO-PHON XL</b> 1 ICU + 1 ICFK + 1 COMBIPHON®</p>	1	5 - 13	<b>21272</b>	
 	<p><b>DURAVENT® MULTI</b> 1 ICU + 1 ICFK + 1 COMBIPHON® O<sub>2</sub></p>	1	5 - 13	<b>21055</b>	
 	<p><b>DURAVENT® SOLO-ADAPTER KOMBI</b></p>	5	5 - 13	<b>11990</b>	
 	<p><b>DURAVENT® SOLO-ADAPTER UNI</b></p>	5	5 - 13	<b>11992</b>	
 	<p><b>DURAVENT® PLUG</b></p>	5	5 - 13	<b>11995</b>	













	Description	Unit / pcs	Size	REF	
 	<b>DURAPUR® LINGO-PHON</b> 1 ICF		7 - 12	<b>23032</b>	
 	<b>DURAPUR® IC</b>		7 - 12	<b>23801</b>	
 	<b>DURAPUR® ICU</b>		7 - 12	<b>23811</b>	
 	<b>DURAPUR® ICK</b>		7 - 12	<b>23881</b>	
 	<b>DURAPUR® ICFU</b>		7 - 12	<b>23816</b>	
 	<b>DURAPUR® ICF</b>		7 - 12	<b>23807</b>	
 	<b>DURAPUR® ICFK</b>		7 - 12	<b>23885</b>	
 	<b>DURAPUR® PHON ICF</b>		7 - 12	<b>23808</b>	
 	<b>DURAPUR® PHON ICF-O2</b>		7 - 12	<b>23809</b>	






	Description	Unit / pcs	Size	REF	
	<b>SENSOTRACH® Foam OVAL</b>	10	oval	<b>30790</b>	
	<b>SENSOTRACH® Foam PRO OVAL</b>	10	oval	<b>30791</b>	
	<b>SENSOTRACH® Foam ROUND</b>	10	round	<b>30792</b>	
	<b>SENSOTRACH® Foam PRO ROUND</b>	10	round	<b>30793</b>	



























# TRACHEOSTOMY TUBES - CUFF

Cuffed tracheostomy tubes are used to seal the cannula against the tracheal wall. This is a requirement for artificial ventilation and/or to prevent accidental aspiration of secretions.























The soft cuff with very thin walls is typically designed as a high volume/low pressure cuff (hv/lp) and fits snugly into the trachea to ensure a reliable seal. The inflation line is integrated into the wall of the (outer) cannula thus providing a rounded and smooth outer shape to avoid friction and irritations of the delicate mucosa, particularly around the tracheostoma.













	Description	Unit / pcs	Size	REF	
 	<b>DURACUFF® VARIO</b>	1	7 - 12	<b>13000</b>	
 	<b>M</b> <b>DURACUFF® VARIO M</b>	1	7 - 12	<b>13100</b>	
 	<b>XL</b> <b>DURACUFF® VARIO XL</b>	1	7 - 12	<b>13200</b>	
 	<b>DURACUFF® UNI-VARIO</b> 1 ICU + 1 ICV	1	7 - 12	<b>13002</b>	
 	<b>M</b> <b>DURACUFF® UNI-VARIO M</b> 1 ICU + 1 ICV	1	7 - 12	<b>13102</b>	
 	<b>XL</b> <b>DURACUFF® UNI-VARIO XL</b> 1 ICU + 1 ICV	1	7 - 12	<b>13202</b>	

	Description	Unit / pcs	Size	REF	
 <p>M</p> <p>      </p>	<b>DURACUFF® UNI M</b> 1 ICU-O <sub>2</sub>	1	7 - 12	<b>13151</b>	
	<b>DURACUFF® CLIP</b>	1		<b>13299</b>	
 <p>     </p>	<b>SPIRAFLEX® UNI CUFF</b>	1	6 - 11	<b>17503</b>	
 <p>      </p>	<b>SPIRAFLEX® MRT UNI CUFF</b>	1	6 - 11	<b>17520</b>	
 <p>      </p>	<b>SPIRAFLEX® SHORT UNI CUFF</b> 2 ICX	1	6 - 11	<b>17507</b>	
 <p>XL</p> <p>      </p>	<b>SPIRAFLEX® UNI CUFF XL</b>	1	7 - 11	<b>17519</b>	
 <p>   </p>	<b>SPIRAFLEX® ICX</b>	3	6 - 11	<b>17543</b>	
 <p>S</p> <p>   </p>	<b>SPIRAFLEX® SHORT ICX</b>	3	6 - 11	<b>17516</b>	
 <p>XL</p> <p>   </p>	<b>SPIRAFLEX® ICX XL</b>	3	7 - 11	<b>17549</b>	

	Description	Unit / pcs	Size	REF	
 <p>Icons:    <b>BE</b> <b>STERILE</b></p>	<b>TRACHEOTEC® VARIO CUFF</b>	1	3 - 10	18920	
 <p>Icons:    <b>BE</b> <b>STERILE</b></p>	<b>TRACHEOTEC® PRO VARIO CUFF</b>	1	6 - 10	18922	
 <p>Icons:  <b>STERILE</b></p>	<b>TRACHEOTEC® PRO IC</b>	3	6 - 10	18933	
 <p>Icons:    <b>BE</b> <b>STERILE</b></p>	<b>TRACHEOTEC® SILC VARIO CUFF</b>	1	5 - 9	18850	
 <p>Icons:    <b>BE</b> <b>STERILE</b></p>	<b>TRACHEOSILC® VARIO CUFF</b>	1	5 - 9,5	18830	
 <p>Icons:    <b>BE</b> <b>STERILE</b></p>	<b>DURATWIX® UNI CUFF</b>	1	7 - 10	19601	
 <p>Icons:    <b>BE</b> <b>STERILE</b></p>	<b>DURATWIX® UNI CUFF short</b>	1	7 - 10	19602	<div style="border: 1px solid black; padding: 2px; display: inline-block;">S</div>






	Description	Unit / pcs	Size	REF	
 	<b>DURATWIX® UNI-VARIO CUFF</b> 1 ICU + 1 ICV	1	6 - 10	<b>19611</b>	
 	<b>DURATWIX® UNI-VARIO CUFF short</b> 1 ICU + 1 ICV	1	7 - 10	<b>19612</b>	
 	<b>DURATWIX® IC</b>	3	7 - 10	<b>19801</b>	
 	<b>DURATWIX® IC short</b>	3	7 - 10	<b>19802</b>	
 	<b>DURATWIX® ICU</b>	3	6 - 10	<b>19811</b>	
 	<b>DURATWIX® ICU short</b>	3	6 - 10	<b>19812</b>	
 	<b>DURATWIX® ICV</b>	3	6 - 10	<b>19821</b>	
 	<b>DURATWIX® ICV short</b>	3	6 - 10	<b>19822</b>	
 	<b>DURATWIX® ICK</b>	3	7 - 10	<b>19831</b>	
 	<b>DURATWIX® ICK short</b>	3	7 - 10	<b>19832</b>	
 	<b>DURATWIX® ICU contrast</b>	3	6 - 10	<b>19803</b>	

























	Description	Unit / pcs	Size	REF	
 	<b>DURATWIX® SOLO-ADAPTER KOMBI</b>	3	22 mm	<b>19808</b>	
 	<b>DURATWIX® SOLO-ADAPTER VARIO</b>	3	15 mm	<b>19807</b>	
 	<b>DURATWIX® SOLO-ADAPTER UNI</b>	3	15 mm	<b>19806</b>	
	<b>MUCOPROTECT® Pressure gauge</b>	1		<b>19500</b>	
	<b>MUCOPROTECT® Connection tube for cuff pressure gauge</b>	1		<b>19506</b>	
	<b>SENSOTRACH® Foam OVAL</b>	10	oval	<b>30790</b>	
	<b>SENSOTRACH® Foam PRO OVAL</b>	10	oval	<b>30791</b>	
	<b>SENSOTRACH® Foam ROUND</b>	10	round	<b>30792</b>	
	<b>SENSOTRACH® Foam PRO ROUND</b>	10	round	<b>30793</b>	































These tracheostomy tubes combine the features of a cuffed tracheostomy tube with a sieve-fenestrated one: Cuffed tracheostomy tubes are used to seal the cannula against the tracheal wall. This is a requirement for artificial ventilation and/or to prevent accidental aspiration of secretions. The soft cuff with very thin walls is typically designed as a high volume/low pressure cuff (hv/lp) and fits snugly into the trachea to ensure a reliable seal. The inflation line is integrated into the wall of the (outer) cannula thus providing a rounded and smooth outer shape to avoid friction and irritations of the delicate mucosa, particularly around the tracheostoma. Sieve-fenestrated tracheostomy tubes have an opening in the upper (convex) wall of the tube shaft facing the larynx, the so-called sieve-fenestration.















While most conventional tracheostomy tubes have one or two simple holes the Fahl tracheostomy tubes have an elaborate sieve consisting of many small holes to reduce micro-aspiration and the ingrowth of granulation tissue. If an inner cannula is present, it has an oval opening corresponding to the sieve in the outer cannula. Sieve-fenestrated cannulas are used by tracheostomized patients with intact larynx who are able to speak, or by laryngectomized patients with voice prosthesis who still need a tracheostomy tube to stabilize the trachea and keep the tracheostoma open. In that case the sieve supports the speaking facility provided by the voice prosthesis.







Sieve-fenestrated tracheostomy tubes are usually used in combination with a speaking valve.



	Description	Unit / pcs	Size	REF	
 <p>BE STERILE</p>	<b>DURATWIX® UNI-VARIO LINGO-CUFF</b> 1 ICV + 1 ICFU	1	6 - 10	19661	
 <p>S</p> <p>BE STERILE</p>	<b>DURATWIX® UNI-VARIO LINGO-CUFF short</b> 1 ICV + 1 ICFU	1	7 - 10	19662	
 <p>BE STERILE</p>	<b>DURATWIX® UNI LINGO-CUFF</b> 1 ICU + 1 ICFU	1	6 - 10	19621	

	Description	Unit / pcs	Size	REF	
 	DURATWIX® IC	3	7 - 10	19801	
 	DURATWIX® IC short	3	7 - 10	19802	
 	DURATWIX® ICU	3	6 - 10	19811	
 	DURATWIX® ICU short	3	6 - 10	19812	
 	DURATWIX® ICV	3	6 - 10	19821	
 	DURATWIX® ICV short	3	6 - 10	19822	
 	DURATWIX® ICK	3	7 - 10	19831	
 	DURATWIX® ICK short	3	7 - 10	19832	
 	DURATWIX® ICU contrast	3	6 - 10	19803	
 	DURATWIX® PHON ICF	2	7 - 10	19841	
 	DURATWIX® PHON ICF short	2	7 - 10	19842	
 	DURATWIX® ICF	3	7 - 10	19851	

	Description	Unit / pcs	Size	REF	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span> </div> <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® ICF short	3	7 - 10	19852	
 <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® ICFU	3	6 - 10	19861	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span> </div> <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® ICFU short	3	6 - 10	19862	
 <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® ICFV	3	6 - 10	19871	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span> </div> <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® ICFV short	3	6 - 10	19872	
 <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® ICFK	3	7 - 10	19881	
 <div style="display: flex; align-items: center; gap: 5px;"> <span style="border: 1px solid black; padding: 2px;">S</span> </div> <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® ICFK short	3	7 - 10	19882	
 <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® SOLO-ADAPTER KOMBI	3	22 mm	19808	
 <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® SOLO-ADAPTER VARIO	3	15 mm	19807	
 <div style="display: flex; align-items: center; gap: 5px;">   </div>	DURATWIX® SOLO-ADAPTER UNI	3	15 mm	19806	

	Description	Unit / pcs	Size	REF	
 	<b>DURATWIX® PLUG</b>	5	15 mm	19805	
 <span data-bbox="443 521 475 560">S</span> 	<b>SPIRAFLEX® SHORT UNI LINGO-PHON CUFF</b> 1 ICX + 1 ICFX + Speaking Valve	1	6 - 11	17509	
 <span data-bbox="443 813 475 851">S</span> 	<b>SPIRAFLEX® SHORT ICX</b>	3	6 - 11	17516	
 <span data-bbox="443 1003 475 1041">S</span> 	<b>SPIRAFLEX® SHORT ICFX</b>	3	6 - 11	17517	
 	<b>DURACUFF® UNI-VARIO LINGO</b> 1 ICV + 1 ICFU	1	7 - 12	13022	
 <span data-bbox="443 1507 475 1545">M</span> 	<b>DURACUFF® UNI-VARIO LINGO M</b> 1 ICV + 1 ICFU	1	7 - 12	13122	
 <span data-bbox="443 1821 475 1859">XL</span> 	<b>DURACUFF® UNI-VARIO LINGO XL</b> 1 ICV + 1 ICFU	1	7 - 12	13222	











	Description	Unit / pcs	Size	REF	
 <p>BE STERILE</p>	<p><b>DURACUFF®</b>  <b>VARIO-KOMBI LINGO</b>                      1 ICV + 1 ICFK</p>	1	7 - 12	13072	
 <p>M</p> <p>BE STERILE</p>	<p><b>DURACUFF®</b>  <b>VARIO-KOMBI LINGO M</b>                      1 ICV + 1 ICFK</p>	1	7 - 12	13172	
 <p>XL</p> <p>BE STERILE</p>	<p><b>DURACUFF®</b>  <b>VARIO-KOMBI LINGO XL</b>                      1 ICV + 1 ICFK</p>	1	7 - 12	13272	
	<p><b>DURACUFF®</b> CLIP</p>	1		13299	
	<p><b>MUCOPROTECT®</b>                      Pressure gauge</p>	1		19500	
	<p><b>MUCOPROTECT®</b>                      Connection tube for cuff pressure gauge</p>	1		19506	












	Description	Unit / pcs	Size	REF	
	<b>SENSOTRACH® Foam OVAL</b>	10	oval	30790	
	<b>SENSOTRACH® Foam PRO OVAL</b>	10	oval	30791	
	<b>SENSOTRACH® Foam ROUND</b>	10	round	30792	
	<b>SENSOTRACH® Foam PRO ROUND</b>	10	round	30793	












Tracheostomy tubes with suction line are cuffed tubes which offer the option of removing secretions which have accumulated above the inflated cuff. This is achieved by a fine tube which is integrated into the (outer) cannula to retain a rounded and smooth outer shape of the cannula to avoid friction and irritations of the delicate mucosa, particularly around the tracheostoma.

The suction line can be connected to any suction unit. It has an orifice right above the cuff to ensure that all secretions can be removed.






	Description	Unit / pcs	Size	REF	
 	<b>DURACUFF® VARIO SUCTION</b>	1	7 - 12	13040	
 	<b>DURACUFF® VARIO SUCTION M</b>	1	7 - 12	13140	
 	<b>DURACUFF® VARIO SUCTION XL</b>	1	7 - 12	13240	
 	<b>TRACHEOTEC® VARIO CUFF-SUCTION</b>	1	5 - 10	18930	
 	<b>DURACUFF® UNI-VARIO SUCTION</b> 1 ICU + 1 ICV	1	7 - 12	13032	

	Description	Unit / pcs	Size	REF	
 <p><b>M</b></p> <p><b>DURACUFF® UNI-VARIO SUCTION M</b> 1 ICU + 1 ICV</p> <p>  </p> <p>STERILE</p>		1	7 - 12	13132	
 <p><b>XL</b></p> <p><b>DURACUFF® UNI-VARIO SUCTION XL</b> 1 ICU + 1 ICV</p> <p>  </p> <p>STERILE</p>		1	7 - 12	13232	
 <p><b>DURACUFF® VARIO-KOMBI LINGO-PHON SUCTION</b> 1 ICV + 1 ICFK + 1 COMBIPHON®</p> <p>  </p> <p>STERILE</p>		1	7 - 12	13042	
 <p><b>M</b></p> <p><b>DURACUFF® VARIO-KOMBI LINGO-PHON SUCTION M</b> 1 ICV + 1 ICFK + 1 COMBIPHON®</p> <p>  </p> <p>STERILE</p>		1	7 - 12	13142	
 <p><b>XL</b></p> <p><b>DURACUFF® VARIO-KOMBI LINGO-PHON SUCTION XL</b> 1 ICV + 1 ICFK + 1 COMBIPHON®</p> <p>  </p> <p>STERILE</p>		1	7 - 12	13242	
 <p><b>DURACUFF® CLIP</b></p>		1		13299	








	Description	Unit / pcs	Size	REF	
 <p>Icons: CURVED, TRACH, CUFF, SUCTION, BE, STERILE</p>	<b>TRACHEOSILC® VARIO CUFF-SUCTION</b>	1	5 - 9,5	18840	
 <p>Icons: CURVED, TRACH, CUFF, SUCTION, BE, STERILE</p>	<b>TRACHEOTEC® PRO VARIO CUFF-SUCTION 2 IC</b>	1	6 - 10	18932	
 <p>Icon: STERILE</p>	<b>TRACHEOTEC® PRO IC</b>	3	6 - 10	18933	
	<b>MUCOPROTECT® Pressure gauge</b>	1		19500	
	<b>MUCOPROTECT® Connection tube for cuff pressure gauge</b>	1		19506	
	<b>SENSOTRACH® Foam OVAL</b>	10	oval	30790	
	<b>SENSOTRACH® Foam PRO OVAL</b>	10	oval	30791	
	<b>SENSOTRACH® Foam ROUND</b>	10	round	30792	
	<b>SENSOTRACH® Foam PRO ROUND</b>	10	round	30793	

The voice of a laryngectomized patient can be restored by placing a voice prosthesis during or after laryngectomy between the esophagus and the trachea. It restores two vital functions of the larynx: It facilitates a continuous air flow from the lungs via the lower trachea, through the voice prosthesis and the esophagus into the oral cavity. The exhaled air can be used for phonation. And it seals the trachea against the esophagus during eating and drinking. Thus no accidental aspiration of food can occur.





There are basically two types of voice prostheses: Indwelling prostheses which need to be exchanged by the physician and non-indwelling prostheses which may be exchanged by the patient himself. A voice prosthesis needs to be replaced for various reasons. E.g. after a certain time the valve may not close properly anymore due to Candida growth or other sediments which may lead to aspiration or an enhanced air flow resistance. For details on this kind of surgical voice rehabilitation a specialized physician should be consulted.

	Description	Unit / pcs	Size	REF *	
	<b>Blom-Singer® CLASSIC Indwelling Voice Prosthesis</b>	1	4 - 20 16 Fr.  4 - 20 20 Fr.	<b>IN16XX -NS</b>  <b>IN20XX -NS</b>	
 <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px;">STERILE</div>	<b>Blom-Singer® CLASSIC Indwelling Voice Prosthesis</b>  <b>STERIL</b>	1	8 - 14 16 Fr.  8 - 14 20 Fr.	<b>IN16XX -SO</b>  <b>IN20XX -SO</b>	
	<b>Blom-Singer® ADVANTAGE Indwelling Voice Prosthesis</b>  <b>Soft Valve Assembly</b>	1	4 - 14 16 Fr.	<b>AS16XX- NS</b>	
	<b>Blom-Singer® ADVANTAGE Indwelling Voice Prosthesis</b>  <b>Hard Valve Assembly</b>	1	4 - 14 20 Fr.	<b>AD20XX- NS</b>	
	<b>Blom-Singer® DUAL VALVE Indwelling Voice Prosthesis</b>	1	6 - 14 20 Fr.	<b>DV20XX- NS</b>	



\* xx = Size voice prosthesis (e.g. IN16XX-NS in size 04 = IN1604-NS)





	Description	Unit / pcs	Size	REF *	
	<b>Blom-Singer® DUAL VALVE Indwelling LF (Large Flange) Voice Prosthesis</b>	1	6 - 14 20 Fr.	<b>DV20XX- LF</b>	
	<b>Blom-Singer® CLASSIC Indwelling SL (Special Length) Voice Prosthesis</b>	1	5 - 9 16 Fr.  5 - 9 20 Fr.	<b>IN16XX- SL  IN20XX- SL</b>	
	<b>Blom-Singer® CLASSIC Indwelling LESL (Large Esophageal Flange/ Special Length) Voice Prosthesis</b>	1	5 + 7 16 Fr.  5 + 7 20 Fr.	<b>IN16XX- LESL  IN20XX- LESL</b>	
	<b>Blom-Singer® CLASSIC Indwelling LEF (Large Esophageal Flange) Voice Prosthesis</b>	1	4 - 14 16 Fr.  4 - 14 20 Fr.	<b>IN16XX -LEF  IN20XX -LEF</b>	
	<b>Blom-Singer® CLASSIC Indwelling LF (Esophageal + Tracheal Flange) Voice Prosthesis</b>	1	4 - 14 16 Fr.  4 - 14 20 Fr.	<b>IN16XX- LF  IN20XX- LF</b>	
	<b>Blom-Singer® CLASSIC Indwelling IR (Increased Resistance) Voice Prosthesis</b>	1	6 - 12 16 Fr.  6 - 12 20 Fr.	<b>IN16XX- IR  IN20XX- IR</b>	
	<b>Blom-Singer® CLASSIC Indwelling LEIR (Large Esophageal Flange / Increased Resistance) Voice Prosthesis</b>	1	4 - 8 16 Fr.  4 - 8 20 Fr.	<b>IN16XX -LEIR  IN20XX -LEIR</b>	

\* xx = Size voice prosthesis (e.g. IN16XX-NS in size 04 = IN1604-NS)

	Description	Unit / pcs	Size	REF *	
	<b>Blom-Singer® CLASSIC Indwelling LETO (Large Esophageal Flange/ TEP Occluder) Voice Prosthesis</b>	1	4 - 14 16 Fr.	<b>IN16XX- LETO</b>	
			4 - 14 20 Fr.	<b>IN20XX- LETO</b>	
	<b>Blom-Singer® DUCKBILL 16 Fr. Voice Prosthesis</b>  <b>patient-changeable</b>	1	6	<b>DB16-006</b>	
			8	<b>DB16-008</b>	
			10	<b>BE6000</b>	
			12	<b>DB16-012</b>	
			14	<b>BE6001</b>	
			18	<b>BE6002</b>	
	<b>Blom-Singer® LOW PRESSURE 16 Fr. Voice Prosthesis</b>  <b>patient-changeable</b>	1	6	<b>BE6009</b>	
			8	<b>LP16-008</b>	
			10	<b>BE6010</b>	
			12	<b>LP16-012</b>	
			14	<b>BE6011</b>	
			18	<b>BE6012</b>	
			22	<b>BE6013</b>	
			25	<b>BE6014</b>	
28	<b>BE6015</b>				
	<b>Blom-Singer® LOW PRESSURE 20 Fr. Voice Prosthesis</b>  <b>patient-changeable</b>	1	6	<b>BE2014</b>	
			8	<b>LP20-008</b>	
			10	<b>BE2018</b>	
			12	<b>LP20-012</b>	
			14	<b>BE2022</b>	
			18	<b>BE2026</b>	
			22	<b>BE2030</b>	
25	<b>BE2033</b>				

\* xx = Size voice prosthesis (e.g. IN16XX-NS in size 04 = IN1604-NS)

	Description	Unit / pcs	Size	REF	
	<p><b>Blom-Singer® Indwelling Plug Insert</b></p> <p>for CLASSIC Indwelling voice prostheses 4 - 25 mm</p>	1	16 Fr. 20 Fr.	IN4069 IN4070	
	<p><b>Blom-Singer® Tracheoesophageal Puncture Dilator</b></p> <p>for 16 + 20 Fr. voice prostheses</p>	1	18 Fr. 22 Fr.	BE6050 BE2050	
	<p><b>Blom-Singer® Dilator-Sizer</b></p> <p>for 16 + 20 Fr. voice prostheses</p>	1	18 Fr. 22 Fr.	BE2051 BE2052	
	<p><b>Blom-Singer® Replacement Gel Caps</b></p>	90 90 60	16 Fr. 20 Fr. 22 Fr.	BE3190 BE3290 BE3295	
	<p><b>LARYVOX® DUO BRUSH incl. Clip</b></p>	6	6 - 18	29110	
	<p><b>LARYVOX® Stoma Light</b></p>	1		69101	
	<p><b>LARYVOX® PLUG</b></p> <p>for all common voice prostheses 20, 22 + 22.5 Fr.</p>	1		25800	

	Description	Unit / pcs	Size	REF	
	<p><b>LARYVOX® INSERTER</b></p> <p>for different voice prostheses</p>	1		25900	
	<p><b>LARYVOX® INSERTER EXPERT</b></p> <p>for 16 Fr. / DUAL VALVE indwelling voice prostheses (20 Fr).</p>	1		25910	
	<p><b>LARYVOX® INSERTER EXPERT SMALL</b></p> <p>for voice prostheses of 4 mm length</p>	1		25915	
	<p><b>LARYVOX® INSERTER SPECIAL</b></p> <p>application with the gel capsule technique</p>	1		25920	



Speaking valves facilitate speaking for tracheotomized patients with a partly or fully intact larynx. They allow the patients to speak without finger occlusion of the tracheostoma. They open during inhalation and close during exhalation thereby guiding the exhaled air through the larynx, the vocal cords and the mouth cavity. Thus the patient is able to speak. They also provide enhanced hygiene since finger occlusion is no longer necessary.

- **Speaking valves with 15 mm standard connector**

The plastic housing is equipped on the patient side with a central opening with a 15 mm internal diameter, thereby allowing connection to a 15 mm standard connector.

- **Speaking valves with 22 mm connector**

In addition, there is also the variant of a 22 mm connector, which is compatible with all commonly used 22 mm receptacle systems.

- **Speaking valves fitted with a multi-adapter** combine both variants (15 and 22 mm).

- **Low resistance speaking valves** are fitted with a laterally suspended valve membrane.

At rest, the valve membrane is slightly opened. Only slight blowing pressure is necessary to close the speaking valve.

- The **HUMIDOPHONE® speaking valve** features an additional filter medium which serves the thermal and humidity exchange of the breathing air and its filtration.

	Description	Unit / pcs	Size	REF	
 	COMBIPHON®	1	15 mm 22 mm	27131	
  	COMBIPHON® O <sub>2</sub>	1	15 mm 22 mm	27132	
	COMBIPHON® Low Resistance	1	15 mm 22 mm	27133	
 	COMBIPHON® O <sub>2</sub> Low Resistance	1	15 mm 22 mm	27134	
	COMBIPHON® SLIM	1	15 mm	27108	

	Description	Unit / pcs	Size	REF	
 <small>O<sub>2</sub></small>	COMBIPHON® SLIM O <sub>2</sub>	1	15 mm	27109	
	HUMIDOPHONE®	30	22 mm	46480	
	HUMIDOPHONE® PLUS	30	15 mm 22 mm	46487	
 <small>O<sub>2</sub></small>	HUMIDOPHONE® PLUS O <sub>2</sub>	30	15 mm 22 mm	46489	

These tracheostoma valves facilitate speaking for laryngectomized patients with a voice prosthesis. They allow the patients to speak without finger occlusion of the tracheostoma. They open during inhalation and close during exhalation thereby guiding the exhaled air through the voice prosthesis, the esophagus and the mouth cavity. Thus even a laryngectomized patient is again able to speak. They also provide enhanced hygiene since finger occlusion is no longer necessary.

Tracheostoma valves come with 22mm adapters for attachment to stoma buttons and base plates.

	Description	Unit / pcs	Size	REF	
	<b>LARYVOX® Hands-Free Valve Low Resistance</b>	1	22 mm	<b>21984</b>	
	<b>LARYVOX® Hands-Free Valve Kombi HME Low Resistance</b>	1	22 mm	<b>21982</b>	
	<b>LARYVOX® Hands-Free Valve Kombi HME Normal Resistance</b>	1	22 mm	<b>21983</b>	
	<b>LARYVOX® Hands-Free Valve TrachiNaze® Normal Resistance</b>	1	22 mm	<b>21981</b>	
	<b>LARYVOX® Hands-Free Valve TrachiNaze® Low Resistance</b>	1	22 mm	<b>21980</b>	
	<b>LARYVOX® Hands-Free Set</b> 1x LARYVOX® Hands-Free Valve Kombi HME Low Resistance 6x LARYVOX® DUO BRUSH - 10 5x LARYVOX® HME 5x LARYVOX® Tape COMFORT 5x OPTIGARD® 5x OPTICLEAR®	1	22 mm	<b>21986</b>	

The HME (Heat-Moisture-Exchanger) filter cassettes condition the air, protect the lower airways, improve the quality of speech and are simple to use. A tracheostomy results in a loss of mouth/nose breathing, this means that the filtering, warming and humidification of the air has to happen in a different way. To prevent irritation of the airways, due to dry air, dust or harmful substances contained in the air, a filter is required. After surgery, tracheostomised/laryngectomised patients frequently suffer from increased production of mucous and coughing. These symptoms can impair speech rehabilitation and breathing. The regular use of HME filter cassettes with Heat and Moisture Exchangers alleviates the symptoms and increases the quality of life.

### Comfortable tracheostoma closure

The housing of the HME filter cassettes facilitates speaking for laryngectomised patients with a voice prosthesis. The tracheostoma is closed for speaking with the simple touch of a finger, even through clothing. As soon as the finger is removed, the air can pass through the HME to facilitate unimpaired inhalation and exhalation.




### Inconspicuous

Thanks to their low profile the HME filter cassettes are very inconspicuous. They are either connected to tracheostomy tubes or to stoma buttons or to special base plates. The user must take care to select the right connector/adaptor. Most tracheostomy tubes have a 15mm connector, stoma buttons and base plates have a 22mm adapter. Fahl HMEs are available with either 15mm connectors or 22mm adapters or the unique combi-adaptor which facilitates the connection to both types.

### Base Plates

are used by laryngectomized patients. They are self-adhesive carrier systems for HMEs or tracheostoma valves (speaking valves for laryngectomized patients). We offer many types of base plates with various shapes and sizes, made of different materials and with different adhesive strengths. They are selected according to the anatomy of the tracheostoma and the type of skin in order to achieve a long wearing period with an optimal attachment.

An HME or tracheostoma valve may be fitted in the central opening with its 22mm adapter. They can be changed separately from the base plate.

	Description	Unit / pcs	Size	REF	
	LARYVOX® HME	30	22 mm	49800	
	LARYVOX® HME HighFlow	30	22 mm	49810	
	LARYVOX® O <sub>2</sub> HME	30	15 mm 22 mm	49802	

O<sub>2</sub>

	Description	Unit / pcs	Size	REF	
	<p><b>LARYVOX® EXTRA HME</b></p>	30	22 mm	<b>49860</b>	
	<p><b>LARYVOX® EXTRA HME HighFlow</b></p>	30	22 mm	<b>49861</b>	
	<p><b>LARYVOX® EXTRA HME Medium</b></p>	30	22 mm	<b>49862</b>	
	<p><b>LARYVOX® EXTRA HME Sport</b></p>	30	22 mm	<b>49863</b>	
	<p><b>LARYVOX® STYLE HME HighFlow SET</b></p> <p>30x LARYVOX® STYLE HME HighFlow            9x LARYVOX® STYLE TOP lamella (3x gold / 3x silver / 3x blue metallic)            9x LARYVOX® STYLE CAP round (3x gold / 3x silver / 3x blue metallic)</p>	1	22 mm	<b>49880-90</b>	

	Description	Unit / pcs	Size	REF	
	<p><b>LARYVOX® STYLE HME SET</b></p> <p>30x LARYVOX® STYLE HME 9x LARYVOX® STYLE TOP lamella (3x gold / 3x silver / 3x blue metallic) 9x LARYVOX® STYLE CAP round (3x gold / 3x silver / 3x blue metallic)</p>	1	22 mm	<b>49881-90</b>	
	<p><b>LARYVOX® STYLE HME HighFlow</b></p> <p>suitable for use with other LARYVOX® STYLE products</p>	30	22 mm	<b>49880</b>	
	<p><b>LARYVOX® STYLE HME</b></p> <p>suitable for use with other LARYVOX® STYLE products</p>	30	22 mm	<b>49881</b>	
	<p><b>LARYVOX® STYLE TOP lamella - Set</b></p> <p>3x gold 3x silver 3x blue metallic</p>	9		<b>49830-90</b>	
	<p><b>LARYVOX® STYLE TOP lamella - blue metallic</b></p>	3		<b>49830-01</b>	
	<p><b>LARYVOX® STYLE TOP lamella - gold</b></p>	3		<b>49830-02</b>	
	<p><b>LARYVOX® STYLE TOP lamella - silver</b></p>	3		<b>49830-03</b>	
	<p><b>LARYVOX® STYLE TOP lamella - beige</b></p>	3		<b>49830-04</b>	







	Description	Unit / pcs	Size	REF	
	<p><b>LARYVOX® STYLE CAP round - Set</b></p> <p>3x gold 3x silver 3x blue metallic</p>	9		49820-90	
	<p><b>LARYVOX® STYLE CAP round - blue metallic</b></p>	3		49820-01	
	<p><b>LARYVOX® STYLE CAP round - gold</b></p>	3		49820-02	
	<p><b>LARYVOX® STYLE CAP round - silver</b></p>	3		49820-03	
	<p><b>LARYVOX® STYLE CAP round - beige</b></p>	3		49820-04	
	<p><b>LARYVOX® TOUCH HME HighFlow - blue metallic</b></p>	30	22 mm	49850-01	
	<p><b>LARYVOX® TOUCH HME HighFlow - gold</b></p>	30	22 mm	49850-02	
	<p><b>LARYVOX® TOUCH HME HighFlow - silver</b></p>	30	22 mm	49850-03	
	<p><b>LARYVOX® TOUCH HME HighFlow - beige</b></p>	30	22 mm	49850-04	


	Description	Unit / pcs	Size	REF	
	<b>LARYVOX® TOUCH HME HighFlow - clear</b>	30	22 mm	<b>49850-13</b>	
	<b>LARYVOX® TOUCH HME HighFlow Set</b>  5x blue metallic 5x gold 5x silver 5x beige 5x clear	25		<b>48830</b>	
	<b>LARYVOX® PRO HME HighFlow</b>	30	22 mm	<b>49840</b>	
	<b>LARYVOX® PRO HME</b>	30	22 mm	<b>49841</b>	
	<b>LARYVOX® TAPE CONVEX</b>	10	oval	<b>48500</b>	
	<b>LARYVOX® TAPE STANDARD</b>	15	round	<b>48100</b>	
	<b>LARYVOX® TAPE STANDARD</b>	15	oval	<b>48200</b>	
	<b>LARYVOX® TAPE STANDARD</b>	15	XL oval	<b>48300</b>	
	<b>LARYVOX® TAPE FLEXIBLE</b>	15	round	<b>48120</b>	



	Description	Unit / pcs	Size	REF	
	LARYVOX® TAPE FLEXIBLE	15	oval	48220	
	LARYVOX® TAPE FLEXIBLE	15	XL oval	48320	
	LARYVOX® TAPE HYDROSOFT	15	round	48130	
	LARYVOX® TAPE HYDROSOFT	15	oval	48230	
	LARYVOX® TAPE HYDROSOFT	15	XL oval	48330	
	LARYVOX® TAPE COMFORT	15	round	48140	
	LARYVOX® TAPE COMFORT	15	oval	48240	
	LARYVOX® TAPE COMFORT	15	XL oval	48340	
	LARYVOX® TAPE HYPOALLERGEN	15	round	48050	
	LARYVOX® TAPE HYPOALLERGEN	15	oval	48150	


	Description	Unit / pcs	Size	REF	
	<p><b>LARYVOX® TAPE HYPOALLERGEN</b></p>	15	XL oval	<b>48250</b>	
	<p><b>LARYVOX® TAPE EXTRA FINE</b></p>	15	recht-eckig	<b>48400</b>	
	<p><b>LARYVOX® Seal</b></p> <p>Cap for use after decannulation. Not suitable for laryngectomized patients.</p>	3	22 mm	<b>48910</b>	
	<p><b>LARYVOX® Connect</b></p>	1	22 > 15 mm	<b>48915</b>	
	<p><b>LARYVOX® Adapt</b></p>	3	15 > 22 mm	<b>48920</b>	
	<p><b>FAHL® Silicone Glue</b></p>	1		<b>29810</b>	
	<p><b>FAHL® Decannulation-Set</b></p> <p>15x LARYVOX® TAPE EXTRA FINE                      1x LARYVOX® Seal                      1x COMBIPHON® O2</p> <p>Not suitable for laryngectomized patients.</p>	1	15 mm 22 mm	<b>48005</b>	

	Description	Unit / pcs	Size	REF	
	<p><b>LARYVOX® TOUCH HME HighFlow Starter Set</b></p> <p>5x 3 LARYVOX® TAPE oval / rectangular                      15x OPTICLEAR®                      15x OPTIGARD®                      15x LARYVOX® TOUCH HME HighFlow - clear                      15x LARYVOX® TOUCH HME HighFlow - beige                      1x LARYVOX® SECUTRACH®</p>	1	22	48820	
	<p><b>LARYVOX® EXTRA HME HighFlow Starter Set</b></p> <p>5x 3 LARYVOX® TAPE oval / rectangular                      15x OPTICLEAR®                      15x OPTIGARD®                      15x LARYVOX® EXTRA HME HighFlow                      1x LARYVOX® SECUTRACH®</p>			48840	
	<p><b>LARYVOX® Best Care For You</b></p> <p>1x ON THE WAY BOX                      1x TAPES-BOX                      1x SHOWER-BOX                      1x NIGHT-BOX                      1x CLEANING-BOX</p>	1		49895	
	<p><b>LARYVOX® Foam S</b></p>	10		48420	
	<p><b>LARYVOX® Foam M</b></p>	10		48430	
	<p><b>LARYVOX® Foam L</b></p>	10		48440	

	Description	Unit / pcs	Size	REF	
 <b>O<sub>2</sub></b>	<b>HUMIDOTRACH®</b>	30	15 mm	<b>46450</b>	
	<b>HUMIDOTWIN®</b>	30	15 mm 22 mm	<b>46460</b>	
 <b>PED</b>	<b>HUMIDOSTOM® Junior</b>	30	15 mm	<b>46850</b>	
 <b>STERILE</b>	<b>HUMIDOFIX®</b>	30	15 mm	<b>46440</b>	
 <b>O<sub>2</sub></b> <b>STERILE</b>	<b>HUMIDUAL®</b>	30	15 mm	<b>46470</b>	
 <b>O<sub>2</sub></b>	<b>HUMIDOSTOM® O<sub>2</sub></b>	30	15 mm	<b>46150</b>	

Tracheal Compresses (Dressings) are basically used for two reasons. They cushion the neck flange of a tracheostomy tube against the skin to avoid lesions, pain and pressure necroses. They further absorb tracheal secretions. Tracheal compresses are available in numerous types, differing in material, number of layers, size, shape and colour. Aluminium coated compresses and foam compresses reduce the adhesion of the compress to dried secretions or clotted blood and thus facilitate an atraumatic replacement.

Compresses with a slit facilitate the replacement without the need to use a pair of scissors to cut them open or to remove the tracheostomy tube. For immediate post-operative use and open wounds a sterile compress should be used.

	Description	Unit / pcs	Size	REF	
	<b>SENSOTRACH® ALU</b>	10	8 x 10 cm	<b>30000</b>	
	<b>SENSOTRACH® ALU SLIT</b>	10	8 x 10 cm	<b>30030</b>	
	<b>SENSOTRACH® ALU SLIT</b>	10	8 x 9 cm	<b>30036</b>	
	<b>SENSOTRACH® UNO</b>	10	9 x 9,8 cm	<b>30520</b>	
	<b>SENSOTRACH® UNO SLIT</b>	10	9 x 9,8 cm	<b>30521</b>	
 <small>PED</small>	<b>SENSOTRACH® UNO PED</b>	10	6 x 7 cm	<b>30522</b>	
 <small>PED</small>	<b>SENSOTRACH® UNO PED SLIT</b>	10	6 x 7 cm	<b>30523</b>	
 <small>PED</small>	<b>SENSOTRACH® UNO PED ALU SLIT</b>	10	6 x 7 cm	<b>30524</b>	

	Description	Unit / pcs	Size	REF	
	<b>SENSOTRACH® DUO</b>	10	9 x 9,8 cm	<b>30608</b>	
	<b>SENSOTRACH® DUO SLIT</b>	10	9 x 9,8 cm	<b>30609</b>	
 <small>PED</small>	<b>SENSOTRACH® DUO PED</b>	10	6 x 7 cm	<b>30622</b>	
 <small>PED</small>	<b>SENSOTRACH® DUO PED SLIT</b>	10	6 x 7 cm	<b>30623</b>	
	<b>SENSOTRACH® DUO ALU</b>	10	9 x 9,8 cm	<b>30662</b>	
	<b>SENSOTRACH® DUO ALU SLIT</b>	10	9 x 9,8 cm	<b>30663</b>	
	<b>SENSOTRACH® DUO SKIN</b>	10	9 x 9,8 cm	<b>30668</b>	
	<b>SENSOTRACH® DUO SKIN SLIT</b>	10	9 x 9,8 cm	<b>30669</b>	
	<b>SENSOTRACH® 3-Plus</b>	10	9 x 9,8 cm	<b>30770</b>	
	<b>SENSOTRACH® 3-Plus SLIT</b>	10	9 x 9,8 cm	<b>30780</b>	










	Description	Unit / pcs	Size	REF	
 <p data-bbox="124 434 199 459">STERILE</p>	SENSOFOAM® - sterile	10	10,0 x 9,0 x 0,5 cm	30860	
 <p data-bbox="124 607 199 631">STERILE</p>	SENSOFOAM® - sterile, small	10	6,5 x 6,5 x 0,5 cm	30861	
 <p data-bbox="124 779 199 804">STERILE</p>	SENSOFOAM® PAD - sterile	10	7,8 x 3,8 cm	30862	
 <p data-bbox="124 949 199 974">STERILE</p>	SENSOFOAM® PAD Professional - sterile	10	7,8 x 3,8 cm	30864	







For reasons of hygiene and to avoid the risk of infection, the tracheostomy tubes should be cleaned thoroughly at least twice a day. And appropriately more often during heavier secretion production. Our range of aids offers a variety of products for gentle and hygienic cleaning, disinfection and preparation of tracheostomy tubes.

Great care should be taken in the handling of tracheostomy tubes, especially when cleaning with a cannula-cleaning brush. Should sharp edges or even cracks ever appear on a tracheostomy tube, on no account continue to use it, as this could lead to injury of the mucous membranes in the trachea. Constant checking avoids danger!

	Description	Unit / pcs	Size	REF	
	OPTIBRUSH® SET	1		31011	
	OPTIBRUSH® CLEAN	1	100 g	31110	
	OPTIBRUSH® CONT	1		31200	
	OPTICIT®	1	250 ml	31180	
	OPTIFLUID® Stoma oil	1	25 ml	31525	
	OPTIFLUID® Stoma oil towel	30		31550	
	OPTIBRUSH®	4	4 - 14	31850	
	OPTIBRUSH® PLUS	4	5 - 14	31855	



	Description	Unit / pcs	Size	REF	
	OPTIBRUSH® Basic	4	6 - 14	31800	
	OPTIBRUSH® Swab	30	small	31910	
	OPTIBRUSH® Swab XL	30	large	31920	
	OPTIFAHL® Stoma cleaning wipes	30	140 x 200 mm	33200	
	OPTIFAHL® Stoma cleaning wipes	30	135 x 200 mm	33210	
	OPTIFAHL® Stoma cleaning wipes box	60		33260	
	OPTICLEAR® Adhesive remover	30		33500	
	OPTIGARD® Protective dressing wipes	30		33600	
	LARYVOX® Distance holder	1		34000	





	Description	Unit / pcs	Size	REF	
	LARYVOX® Tracheospreizer	1		35500	
 <small>PED</small>	LARYVOX® Tracheospreizer PED	1		35501	
	LARYVOX® Tweezers	1		35000	
	LARYVOX® SPECULUM	1		35100	
	FAHL® Gauze Dressings 8-fold	100	10 x 10 cm	30200	
	FAHL® Non-Woven Compress / 4-fold	100	10 x 10 cm	30300	











## TUBE HOLDERS

## AIRWAY MANAGEMENT

We offer an extensive range of tracheostomy tube accessories. A tube holder is indispensable for securely fixing the tube in the tracheostoma. Tube holders are available in a large selection and should be matched to the tracheostomy tube. We provide tube holders, in differing material qualities and fastening techniques.

- H** - with plastic hook
- HM** - with metal hook
- K** - with hook-and-loop fastener

	Description	Unit / pcs	Size	REF	
	<b>KACLIP® HM</b> beige	1	30 - 50 cm	<b>32000</b>	
	<b>KACLIP® H</b> beige	1	30 - 50 cm	<b>32005</b>	
	<b>KACLIP® H</b> white	1	30 - 50 cm	<b>32005-02</b>	
	<b>KACLIP® Clear H</b> transparent	1	29 - 52 cm	<b>32015</b>	
 <small>PED</small>	<b>SUPRAPHIX® 11 H PED</b>	1	15 - 26 cm	<b>32308</b>	
 <small>PED</small>	<b>SUPRAPHIX® 11 K PED</b>	1	18 - 25 cm	<b>32509</b>	
 <small>PED</small>	<b>SUPRAPHIX® 13 H PED</b>	1	17 - 28 cm	<b>32310</b>	
 <small>PED</small>	<b>SUPRAPHIX® 13 K PED</b>	1	20 - 27 cm	<b>32510</b>	
 <small>PED</small>	<b>SUPRAPHIX® 15 H PED</b>	1	19 - 30 cm	<b>32311</b>	

	Description	Unit / pcs	Size	REF	
 <small>PED</small>	SUPRAFIX® 15 K PED	1	22 - 29 cm	32512	
	SUPRAFIX® 18 H	1	22 - 41 cm	32301	
	SUPRAFIX® 18 K	1	25 - 40 cm	32502	
	SUPRAFIX® 25 H	1	35 - 42 cm	32300	
	SUPRAFIX® 25 K	1	32 - 48 cm	32500	
	SUPRAFIX® 30 H	1	34 - 54 cm	32305	
	SUPRAFIX® 30 K	1	37 - 53 cm	32501	
	SUPRAFIX® 36 H	1	40 - 58 cm	32303	
	SUPRAFIX® 36 K	1	43 - 57 cm	32503	
	OPTIFLAUSCH® K	1	24 - 45 cm	32550	

	Description	Unit / pcs	Size	REF	
	OPTIFLAUSCH® K beige	1	24 - 45 cm	32553	
	OPTIFLAUSCH® H	1	24 - 47 cm	32551	
 <small>PED</small>	OPTIFLAUSCH® K PED	1	16 - 25 cm	32555	
 <small>PED</small>	OPTIFLAUSCH® K Junior PED	1	16 - 25 cm	32556	
 <small>PED</small>	OPTIFLAUSCH® Slim K	1	24 - 42 cm	32552	
 <small>PED</small>	OPTIFLAUSCH® Slim K PED	1	16 - 30 cm	32557	
 <small>PED</small>	OPTIFLAUSCH® Slim K Junior PED	1	16 - 30 cm	32558	
	NECKFIX® K	1	45 cm	32620	
	NECKFIX® Blue K	1	45 cm	32610	
 <small>PED</small>	NECKFIX® Blue PED K	1	34 cm	32612	

	Description	Unit / pcs	Size	REF	
	NECKFIX® Beige K PED	1	34 cm	32621	
	NECKFIX® Beige K	1	46 cm	32622	
	NECKFIX® Beige K XL	1	50 cm	32623	
	NECKFIX® Duo K	1	40 - 53 cm	32630	
	NECKFIX® Duo K beige	1	40 - 53 cm	32631	
	NECKFIX® DUO K Plus	50	16 - 45 cm	32629	
	NECKFIX® ENDO Endotracheal tube fixation	10	85 cm	32615	

The tracheostoma reacts particularly sensitive to all outside influence. Since after a tracheostomy or a laryngectomy breathing occurs no longer through the mouth and the nose but through the tracheostoma some of the body's natural protective measures are lost. It is therefore important to attach specific care to the protection of the tracheostoma.

This refers particularly to the humidification, warming and filtering of the inhaled air. The lack of these protective mechanisms may have an adverse effect on the health. For example, a chilled trachea often results in damaging the delicate mucous membranes. These, in turn, perform such important tasks as the transportation of fluids and the movement of foreign bodies out of the upper and lower airways were they normally can be coughed up. The dissication of the tracheal mucosa will cause irritations and increased secretion, partially with the formation of encrustations.





Our products for trachestoma protection replace the necessary protective mechanisms. On the one hand they protect the tracheostoma and trachea from the cold and from drying out – on the other hand they prevent the ingress of dirt particles and absorb leaking secretions.

Without tracheostoma protection the inhaled air would enter through the tracheostoma too cold and too fast. Specific protection items provide a barrier, reducing the speed or the airflow, warming and humidifying the air. At the same time tracheostoma protection products provide visual protection and with it, a certain degree of security.

Our tracheostoma protection products combine aesthetics with comfort. The neckerchiefs and scarves, for example, are offered in fashionable, attractive design patterns.











**Please refer to our colour sample card for the full range.**  
**You can find the current version at [www.fahl.de](http://www.fahl.de).**

- 3L - 3 layers**
- 5L - 5 layers**
- 8L - 8 layers**
- 12L - 12 layers**
- 16L - 16 layers**










	Description	Unit / pcs	Size	REF	
	TRACHEOTEX® BIB 3L	1	34 - 46 cm	40100	
	TRACHEOTEX® BIB 5L	1	34 - 46 cm	40200	
	TRACHEOTEX® BIB 8L	1	34 - 46 cm	40300	
	TRACHEOTEX® BIB 12L	1	34 - 46 cm	40400	

	Description	Unit / pcs	Size	REF	
	TRACHEOTEX® BIB 16L	2	34 - 46 cm	41600-D	
 <small>PED</small>	TRACHEOTEX® BIB 3L PED	1	24 - 33 cm	40110	
 <small>PED</small>	TRACHEOTEX® BIB 5L PED	1	24 - 33 cm	40210	
 <small>PED</small>	TRACHEOTEX® Scarf PED	1	22 - 29 cm	42010	
	BUCHANAN BIB L	10	max. 95 cm 21 x 19 cm	41700	
	BUCHANAN BIB S	10	max. 65 cm 15 x 12 cm	41800	
	TRACHEOTEX® Scarf	1	41 - 47 cm	42000	
	TRACHEOTEX® EASY	1	max. 55 cm	42002	
	TRACHEOTEX® Shirt	1	min. 20 cm	43301	



	Description	Unit / pcs	Size	REF	
	TRACHEOTEX® ROLLI ZIP	1	min. 18 cm	43302-1	
	TRACHEOTEX® ROLLI ZIP FINE	1	min. 18 cm	43302-0	
	TRACHEOTEX® ROLLI KLETT	1	26 - 42 cm	43000-1	
	TRACHEOTEX® ROLLI KLETT FINE	1	26 - 42 cm	43000-0	
	LARYNGOFIX® S beige	10	4 x 5,5 cm	45900	
	LARYNGOFIX® L beige	10	4,8 x 6,7 cm	45950	
	LARYNGOFIX® S white	10	4 x 5,5 cm	45905	
	LARYNGOFIX® L white	10	4,8 x 6,7 cm	45955	
	LARYNGOFIX® PLUS S beige	10	4 x 5,5 cm	45902	
	LARYNGOFIX® PLUS L beige	10	4,8 x 6,7 cm	45952	


	Description	Unit / pcs	Size	REF	
	LARYNGOFIX® PRO S	10	4,0 x 5,5 cm 5 mm	45901	
	LARYNGOFIX® PRO L	10	4,8 x 6,7 cm 5 mm	45956	
	TRACHEOSTOMA-FIX® S beige	10	6 x 6,5 cm	45000	
	TRACHEOSTOMA-FIX® L beige	10	7 x 7,5 cm	45100	
	TRACHEOFIX® S white	10	5,5 x 6 cm	45200	
	TRACHEOFIX® L white	10	7 x 7 cm	45300	
	TRACHEOFIX® S beige	10	5,5 x 6 cm	45400	
	TRACHEOFIX® L beige	10	7 x 7 cm	45500	
	TRACHEOFIX® PLUS S beige	10	5,5 x 6 cm	45505	
	TRACHEOFIX® PLUS L beige	10	7 x 7 cm	45510	

	Description	Unit / pcs	Size	REF	
	LARYNGOFIX® TOUCH L	10	4,8 x 6,7 cm	45911	
	LARYTAPE® S	10	4 x 5,5 cm	45960	
	LARYTAPE® L	10	4,8 x 6,7 cm	45962	
	LARYTAPE® TOUCH S	10	4 x 5,5 cm	45970	
	LARYTAPE® TOUCH L	10	4,8 x 6,7 cm	45972	
	LARYTAPE® XL	10	6,0 x 8,5 cm	45963	
	LARYTAPE® TOUCH XL	10	6,0 x 8,5 cm	45973	
	LARYVOX® SECUTRACH® Shower guard	1	15 mm 22 mm	47400	
	SECUTRACH® Shower guard	1		47000	

Following laryngectomy or tracheotomy, the patient breathes via a tracheostoma (surgical opening of the trachea). When breathing through the nose, the inhaled air is kept moist in a natural manner. After the loss of continuous physiological humidification of the mucous membranes a different method is required in order to prevent severe irritation of the tissue.







In most cases this type of irritation leads to increased secretion, coughing and encrustation. To avoid these consequences, it is essential for the patient to inhale regularly and ensure that his environment are kept at a constant air humidity of approx. 55-60 %.









	Description	Unit / pcs	Size	REF	
	NEBUFIRST®	1		50000	
	NEBUSTEAM®	1		50500	
	NEBUPORT®	1		50100	
	LARYVOX® Nebulizer	1		58000	
	NEBUJUNIOR® Bear	1		50300	
	NEBUJUNIOR® Penguin	1		50400	

	Description	Unit / pcs	Size	REF	
 <div data-bbox="124 488 201 517" style="border: 1px solid black; padding: 2px; width: fit-content;">STERILE</div>	HUMIDOBAC® / HME	50		46830	
	HUMIDOBAC®	50		46831	
	FAHL® Catheter mount XL	1	40 cm	19100	
	FAHL® Catheter mount	1	15 cm	19110	
	FAHL® Nebulizer Set	1		51202	
	FAHL® Oxytube	1		75275	

Tracheal suction devices are important aids for providing care after laryngectomy or tracheotomy. After surgery, patients generally experience increased secretion. This places a great strain on the respiratory system, leading to irritation and coughing. A tracheal suction device is an important aid, helping to remove secretions from the respiratory tract.

The tracheal suction device is also a great help in other situations. In such cases, endotracheal suctioning with the tracheal suction device alleviates breathing, as deep-seated bronchial secretion can be removed professionally. Suction takes place through sterile tracheostoma catheters, which are available in various sizes and types.



	Description	Unit / pcs	Size	REF	
	<b>TRACHEOPORT® COMPACT</b> Battery	1	24 l	<b>63600</b>	
	<b>TRACHEOFIRST® COMPACT</b>	1	24 l	<b>67600</b>	
	<b>TRACHEOPORT® JUNIOR</b>	1	24 l	<b>63800</b>	
	<b>Stoma catheter Metric L</b>	30	50 cm CH 6 - 18	<b>68000</b>	
	<b>Tracheostoma catheter Metric S</b>	30	25 cm CH 10 - 16	<b>68100</b>	
	<b>TRACHFLOW® Metric S</b>	30	25 cm CH 8 - 16	<b>68630</b>	

	Description	Unit / pcs	Size	REF	
 <p data-bbox="121 450 201 472">STERILE</p>	<p data-bbox="521 320 730 387"><b>TRACHFLOW® Metric L</b></p>	<p data-bbox="914 320 954 342">30</p>	<p data-bbox="1078 320 1166 454">50 cm CH 8 - 16</p>	<p data-bbox="1206 320 1286 342"><b>68730</b></p>	
 <p data-bbox="121 680 201 703">STERILE</p>	<p data-bbox="521 528 730 595"><b>TRACHFLOW® SHUNT</b></p>	<p data-bbox="914 528 954 551">50</p>	<p data-bbox="1078 528 1166 663">10 cm CH 8, 10</p>	<p data-bbox="1206 528 1286 551"><b>68800</b></p>	
 <p data-bbox="121 911 201 934">STERILE</p>	<p data-bbox="521 736 794 848"><b>RESPITRACH® Permanent suction catheter</b></p>	<p data-bbox="914 736 930 759">8</p>	<p data-bbox="1078 736 1118 804">10 - 16</p>	<p data-bbox="1206 736 1286 759"><b>68400</b></p>	
 <p data-bbox="121 1120 201 1142">STERILE</p>	<p data-bbox="521 967 882 1001"><b>TRACHFLOW® Finger-Tip</b></p>	<p data-bbox="914 967 954 990">10</p>		<p data-bbox="1206 967 1286 990"><b>60700</b></p>	
	<p data-bbox="521 1176 730 1243"><b>TRACHFLOW® Stepconnector</b></p>	<p data-bbox="914 1176 930 1198">1</p>		<p data-bbox="1206 1176 1318 1198"><b>62000-3</b></p>	
	<p data-bbox="521 1326 858 1359"><b>TRACHFLOW® Line Pro</b></p>	<p data-bbox="914 1326 930 1348">1</p>	<p data-bbox="1078 1326 1166 1348">1,30 m</p>	<p data-bbox="1206 1326 1286 1348"><b>60506</b></p>	
	<p data-bbox="521 1489 850 1556"><b>TRACHFLOW® Silicone suction tube S</b></p>	<p data-bbox="914 1489 930 1512">1</p>	<p data-bbox="1078 1489 1150 1512">28 cm</p>	<p data-bbox="1206 1489 1286 1512"><b>60510</b></p>	
	<p data-bbox="521 1653 799 1720"><b>TRACHFLOW® Line Connect</b></p>	<p data-bbox="914 1653 930 1675">1</p>	<p data-bbox="1078 1653 1150 1675">30 cm</p>	<p data-bbox="1206 1653 1286 1675"><b>60517</b></p>	

After successful laryngectomy, it is most important for the patient to be able to speak as quickly as possible again and to be well understood. That is why we recommend speech therapy by a certified speech therapist.

The speech therapist will usually start speech rehabilitation with the esophageal voice. Parallel to this, speaking with an electronic speech aid can also be a part of the voice rehabilitation. The electronic speech aid generates sound oscillations, these are transmitted to the mouth/pharynx area when the device is slightly pressed to the neck, and can thereby transform the normal speech motions into comprehensible speech.


The VOICETEC® voice amplifier enables excellent amplification of a quiet oesophageal voice and is also well suited for patients/users with vocal chord damage. It can be used whenever a louder conversation volume is desired, to, e.g., be understood even with high levels of background noise. With the help of the VOICETEC® voice amplifier, it is possible, even without strenuous attempts, to speak loudly. The wireless, radio-controlled transmitter unit enables the greatest of mobility. This extends the patient’s comfort and field of application.




	Description	Unit / pcs	Size	REF	
	VOICETEC®	1		77350	
	ElectroLarynx Speech Aid	1		EL1000	



Patients with a tracheostomy unit are generally more prone to emergency situations. On the one hand, the limited or lacking option of speaking can possibly impede effectively calling for help. On the other hand, modified respiration via the tracheostoma cannot be recognised immediately by paramedics or be assigned to appropriate measures.

In order to provide those affected with the greatest possible security in emergencies and the fastest possible emergency care (e.g. air supply), we offer special aids which could prove life-saving in individual cases.

	Description	Unit / pcs	Size	REF	
	<b>LARYVOX® SOS safe-mask</b>	1		75000	
	<b>LARYVOX® SOS mask</b>	1		75010	
	<b>LARYVOX® Smell</b>	1	22 mm	75130	
	<b>LARYVOX® EMERGENCY-SET Laryngectomy</b> 1x LARYVOX® SOS safe-mask 1x LARYVOX® Tweezers 1x LARYVOX® Connect 5x LARYVOX® O <sub>2</sub> Connect 5x LARYVOX® O <sub>2</sub> HME 5x LARYVOX® TAPE COMFORT / XL oval 1x Oxygen connecting tube	1	15 mm 22 mm	48006	
	<b>LARYVOX® EMERGENCY-SET Laryngectomy PRO</b> 1x LARYVOX® SOS safe-mask 1x TRACHEOTEC® VARIO CUFF, Size 7 1x Single-use syringe / 10ml 1x LARYVOX® Tracheospreizer 1x LARYVOX® Connect 5x LARYVOX® O <sub>2</sub> Connect 5x LARYVOX® TAPE COMFORT / XL oval 3x OPTIFLUID® Stoma oil towel 1x Oxygen connecting tube	1	15 mm 22 mm	48007	



	Description	Unit / pcs	Size	REF	
	<p><b>LARYVOX® EMERGENCY-SET Tracheostomy</b></p> <ul style="list-style-type: none"> <li>1x LARYVOX® Tracheospreizer</li> <li>1x HUMIDOTRACH®</li> <li>1x Oxygen connecting tube</li> <li>1x OPTIFLAUSCH® K</li> <li>1x 2-KAM® SLIT, sterile</li> <li>3x OPTIFLUID® Stoma oil towel</li> <li>3x OPTIFAHL® Stoma cleaning wipes</li> <li>1x LARYVOX® SOS safe-mask</li> <li>1x Single-use syringe / 10ml</li> </ul>	1	15 mm	48008	
	<p><b>FAHL® KONIOTOMIE-SET</b></p>	1		75001	
	<p><b>LARYVOX® Alarm</b></p>	1		76000	

The LARCHEL® water therapy device is a proven product, which enables laryngectomised patients to take part in exercise therapy in water and therapeutic swimming. In detail, LARCHEL® offers the following areas of use: underwater massage, aquatic therapy, showering, full bath, smelling and swimming.

The LARCHEL® consists of a flexible corrugated tube, which is firmly combined with a special, blockable tracheostomy tube at one end. A mouthpiece with bite block is located at the other end. This mouthpiece is placed in the user's mouth and firmly surrounded by the lips.

The LARCHEL® is available in different sizes (tracheostomy tube inner diameter). The exact cannula diameter must be selected in order to achieve a good seal during inflation of the tracheostomy tube. The treating physician and an authorised specialist, the aquatic therapy officer of the Association of Laryngectomees, should be consulted as a part of the introduction!

A snorkel is available as addition to the LARCHEL® water therapy device. This is particularly suited for users who experience problems when breathing through their nose.

	Description	Unit / pcs	Size	REF	
	<b>LARCHEL® water therapy device</b>	1	5 - 9,5	79160	
	<b>Snorkel for LARCHEL®</b>	1		79161	

### Blom-Singer® Voice Prosthesis Placement Surgical Kit







The Blom-Singer® Surgical Kit is a complete surgery set for single use and serves the intraoperative primary or secondary puncture for placing a Blom-Singer® Classic™ indwelling voice prosthesis.

Next to the separately packaged voice prosthesis it contains the following components:

- ❶ Catheter (disposable item)
- ❷ Trocar (disposable instrument)
- ❸ Pharynx protector (disposable instrument)


### Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis


The individually adjustable bi-flanged fistula prosthesis was developed for the healthcare management of hypopharyngeal fistulas - in particular, pharyngo-cutaneous fistulas. It offers the user a further innovative option as part of the general professional management of fistulas.

	Description	Unit / pcs	Size	REF	
	<b>Blom-Singer® Voice Prosthesis Placement Surgical Kit</b>	1	16 Fr.	<b>TP1000</b>	
			20 Fr.	<b>TP1001</b>	
	<b>Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis 25 mm</b>  <b>SMALL</b>	1	25 mm	<b>FM1003</b>	
	<b>Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis 38 mm</b>  <b>MEDIUM</b>	1	38 mm	<b>FM1004</b>	
	<b>Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis 50 mm</b>  <b>LARGE</b>	1	50 mm	<b>FM1005</b>	
	<b>Blom-Singer® Fistula Prosthesis Insertion and Removal Catheter</b>  For size SMALL (25 mm) and MEDIUM (38 mm)	1		<b>FM1010</b>	
	<b>FAHL® KONIOTOMIE-SET</b>	1		<b>75001</b>	























**FAHL IS THERE FOR YOU THROUGHOUT EUROPE! - HERE WE ARE REPRESENTED**

---



## SIZING TRACHEOSTOMY TUBES

Tracheostomy tubes are available in different materials and different versions. A good fit of the cannula is essentially determined by two factors, the cannula diameter and the cannula length.

LARYNGOTEC® / PHONTRACH®					DURAPUR® standard / XL					DURAVENT® standard / XL					LARYNGOTEC® Button standard / short				SILVERVENT®									
Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)	Size	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)	
																		0	4.5	6.1	2.5	2.8	50					
																		1	5.0	6.4	3.0	3.2	50					
																		2	5.7	6.8	3.7	3.6	50					
																		3	6.4	7.5	4.4	4.2	55					
																		3.5	4.8	6.2	3.6	2.7	55 / -					
																		4	5.8	7.2	4.4	3.3	55 / -					
																		5	6.6	8.4	5.4	4.3	55 / 95					
																		6	7.8	9.4	6.2	4.9	60 / 95					
																		7	8.8	10.2	7.0	5.5	65 / 100					
																		8	9.8	11.2	8.0	6.5	70 / 105					
																		9	11.2	12.4	9.0	7.5	75 / 105					
																		10	12.2	13.4	10.0	8.5	85 / 110					
																		11	13.2	14.4	11.0	9.5	87 / 115					
																		12	14.2	15.4	12.0	10.2	90 / 120					
																		13	15.2	16.6	13.0	11.2	90 / 125					
																		14	16.6	18.0	14.0	12.2	90 / 125					

## CIRCULAR TEMPLATES

To determine the outer diameter at the cannula tip in mm. From 3-15.5 mm in 0.5 mm steps.



In the selection of tracheostomy tubes of different materials, it should be noted that the material strength could also influence the inner lumen of the cannula. For example, because of the thick-walled material, a non-metal cannula has a smaller inner volume than the thin-walled silver cannula.

# SIZING TRACHEOSTOMY TUBES

As the sizing of the individual cannula manufacturers is based on different measuring principles, these cannot always be directly compared. Exact measurement is often necessary, whereby certain aids can be helpful in measuring, e.g., a flexible measuring tape.

DURACUFF® standard / medium / XL						DURATWIX® standard / short					DURATWIX® JUNIOR PED				SPIRAFLEX® / SPIRAFLEX® MRT					SPIRAFLEX® XL				SPIRAFLEX® short									
Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)	Size	O.D. Tip (mm)	I.D. OC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)		
												3	4.5	3.0	36																		
												3.5	5.3	3.5	39																		
												4	6.0	4.0	42																		
												4.5	6.5	4.5	45																		
												5	7.0	5.0	48																		
												5.5	7.7	5.5	52																		
						6	10.3	11.9	8.0	6.3	82 / -	6	8.3	6.0	55	6	8.4	8.4	6.0	max. 52					6	9.4	9.4	6.0	4.5	61			
7	8.8	10.2	7.0	5.5	65 / 75 / 100	7	10.8	11.9	8.6	7.0	82 / 72					7	10.2	10.2	7.0	max. 68	7	10.0	10.0	7.0	max. 155	7	10.4	10.4	7.0	5.5	67		
8	9.8	11.2	8.0	6.5	70 / 80 / 105	8	11.4	12.5	9.5	8.0	85 / 75					8	11.2	11.2	8.0	max. 88	8	11.0	11.0	8.0	max. 155	8	11.4	11.4	8.0	6.5	77		
9	11.2	12.4	9.0	7.5	75 / 85 / 105	9	12.5	13.7	10.6	9.0	88 / 78					9	12.2	12.2	9.0	max. 117	9	12.0	12.0	9.0	max. 155	9	12.4	12.4	9.0	7.5	77		
10	12.2	13.4	10.0	8.5	85 / 95 / 110	10	13.8	15.0	11.6	10.0	91 / 81					10	13.2	13.2	10.0	max. 118	10	13.0	13.0	10.0	max. 155	10	13.4	13.4	10.0	8.5	88		
11	13.2	14.4	11.0	9.5	87 / 97 / 115											11	14.2	14.2	11.0	max. 123	11	14.0	14.0	11.0	max. 155	11	14.4	14.4	11.0	9.5	88		
12	14.2	15.4	12.0	10.2	90 / 100 / 120																												

O.D. neck flange = Outer diameter behind the neck flange · O.D. Tip = Outer diameter of the outer cannula at the tip of the cannula · I.D. = Inner diameter at the tip of the cannula · O-cannula = Outer cannula · I-cannula = Inner cannula · Length = Length over the outer curve

## CIRCULAR TEMPLATES

To determine the outer diameter at the cannula tip in mm. From 3 - 15.5 mm in 0.5 mm steps.



In comparison, the inner diameter at the cannula tip will therefore play a role. In some tracheostomy tubes, the lengths of the middle radius or the middle length are also stated. These are the mean bending radius of the cannula tube, i.e., the radius of the (geometric) middle layer of the tube bend.

# SIZING TRACHEOSTOMY TUBES

The following overview shows how the cannula diameter and length can be determined. It also illustrates which measurements are most commonly used.

TRACHEOTEC® TRACHEOTEC® PRO						TRACHEOTEC® SILC				TRACHEOSILC®				FAHL® BIESALSKI								
Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	Length (mm)	Length (Midline) (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	Length (mm)	Size	O.D. Tip (mm)	O.D. Shield (mm)	I.D. OC (mm)	I.D. IC (mm)	Length (mm)	
3	4.7	4.7	3.0	47.2	43.7																	
3.5	5.3	5.3	3.5	50.6	46.7																	
4	6.0	6.0	4.0	57.1	52.7											4	5.6	8.4	4.0	2.3	44	
4.5	6.7	6.7	4.5	57.7	52.7																	
5	7.3	7.3	5.0	61.1	55.7	5	7.3	7.3	5.0	59.0	5	7.3	7.3	5.0	57.0	5	6.5	9.0	5.0	3.3	46	
5.5	8.0	8.0	5.5	63.6	57.7																	
6	8.7	8.7	6.0	67.1	60.7	6	8.7	8.7	6.0	70.0	6	8.7	8.7	6.0	63.0	6	7.7	9.7	6.0	3.9	52	
6.5	9.3	9.3	6.5	73.5	66.6																	
7	10.0	10.0	7.0	77.0	69.6	7	10.0	10.0	7.0	75.0	7	10.0	10.0	7.0	71.0	7	8.6	11.0	7.0	5.0	60	
7.5	10.7	10.7	7.5	81.5	73.6						7.5	10.7	10.7	7.5	73.0							
8	11.3	11.3	8.0	87.9	79.5	8	11.0	11.0	8.0	81.0	8	11.0	11.0	8.0	75.0	8	9.7	12.0	8.0	5.9	66	
8.5	12.0	12.0	8.5	90.4	81.5						8.5	11.7	11.7	8.5	78.0							
9	12.7	12.7	9.0	94.9	85.5	9	12.3	12.3	9.0	92.0	9	12.3	12.3	9.0	80.0	9	11.5	14.0	9.0	7.0	71	
											9.5	13.3	13.3	9.5	83.0							
10	14.0	14.0	10.0	105.3	95.0											10	12.1	15.0	10.0	7.9	75	
																11	12.4	16.0	11.0	9.2	78.5	
																12	13.9	17.0	12.0	9.8	77	

O.D. neck flange = Outer diameter behind the neck flange · O.D. Tip = Outer diameter of the outer cannula at the tip of the cannula · I.D. = Inner diameter at the tip of the cannula · O-cannula = Outer cannula · I-cannula = Inner cannula  
Length = Length over the outer curve

## CIRCULAR TEMPLATES

To determine the outer diameter at the cannula tip in mm. From 3 - 15.5 mm in 0.5 mm steps.



13



13.5



14



14.5



15



15.5

The table is intended to facilitate the determination of the outer diameter at the cannula tip. The diameter can simply be determined by placing the cannula (tip) on the template.



# PICTOGRAM LEGEND

Pictogram	Description
<b>Outer cannula</b>	
	Outer cannula with no inner cannula
	Outer cannula with one inner cannula
	Outer cannula with two inner cannulas
	Outer cannula with three inner cannulas
	Outer cannula including 15 mm-connector (UNI)
	Outer cannula including 22 mm-connector (KOMBI)
	Outer cannula including 15 mm-swivel connector (VARIO)
	Stoma button
	Decannulation plug
	Several small fenestration holes in a specific arrangement in the center of the tube. is referred to as „sieve“ (LINGO)
	Outer cannula fenestrated
	Outer cannula with a slit
	With low-pressure cuff (CUFF)
	With suction device ( SUCTION)
	Speaking valve (PHON)
<b>Inner cannula</b>	
	Inner cannula with low profile
	Inner cannula with 15 mm-connector (UNI)
	Inner cannula with 15 mm-swivel connector (VARIO)
	Inner cannula with 22 mm-connector (KOMBI)
	Inner cannula fenestrated
	Inner cannula flexible
	With oxygen support
<b>General</b>	
	Length short
	Length medium
	Length XL
	MRI suitable
	Specially for children
	For artificial ventilation
	Sterile

## LIST OF ABBREVIATIONS

Abbreviation	Description
LINGO	Sieve-fenestrated
CUT	With a slit
CUFF	With low-pressure cuff
SUCTION	With suction device
MULTI	Multi function
XL	Length XL
M	Length medium
Short	Length short
PHON	With speaking valve
UNI	With 15 mm-connector
VARIO	With 15 mm-swivel connector
KOMBI	With 22 mm-connector
FIX	Button with fastening eyelets
SV	Speaking valve
HP	HUMIDOPHONE®
CP	COMBIPHON®

Abbreviation	Description
SZ	Size
IC	Inner cannula
ICF	Inner cannula. fenestrated
ICU	IC with 15 mm-connector (UNI)
ICFU	IC with 15mm-connector (UNI). fenestrated
ICV	IC with 15 mm-swivel connector (VARIO)
ICFV	IC with 15 mm-swivel connector (VARIO). fenestrated
ICK	IC with 22 mm-connector (KOMBI)
ICFK	IC with 22 mm-connector (KOMBI). fenestrated
ICX	IC flexible
ICFX	IC flexible. fenestrated
O2	IC with oxygen support

FAHL

**Andreas Fahl  
Medizintechnik-Vertrieb GmbH**  
August-Horch-Str. 4a  
51149 Cologne - Germany  
Phone +49 (0) 22 03/29 80-0  
Fax +49 (0) 22 03/29 80-100  
mail [info@fahl.de](mailto:info@fahl.de)  
[www.fahl.de](http://www.fahl.de)

**Local branch Berlin**  
Am Borsigturm 62  
13507 Berlin - Germany  
Phone +49 (0) 30/43 73 139-0  
Fax +49 (0) 30/43 73 139-35

**Local branch Munich**  
Aurbacher Straße 2  
81541 Munich - Germany  
Phone +49 (0) 89/58 90 90 9-0  
Fax +49 (0) 89/58 90 90 9-20

**Overseas branch Austria**  
Triester Straße 14 / Bürohaus 2  
2351 Wiener Neudorf - Austria  
Phone +43 (0) 22 36/89 28 38-0  
Fax +43 (0) 22 36/89 28 38-14

**Overseas branch Switzerland**  
Cherstrasse 4  
8152 Glattbrugg - Switzerland  
Phone +41 (0) 448293060  
Fax +41 (0) 448293061



[www.fahl.de](http://www.fahl.de)

\* Andreas Fahl Medizintechnik-Vertrieb GmbH is exclusive sales partner for Blom-Singer® products in some countries in Europe.  
For more information see our web site [www.fahl.de](http://www.fahl.de).